

[Cookie policy](#)

This website uses cookies to ensure you get the best experience on our website. [Cookie policy](#)

[Privacy policy](#)

[Decline](#) [Allow cookies](#)

Tools, Technologies and Training for Healthcare Laboratories



Westgard QC



[My Cart](#) | [Check Out](#) | [Login](#)

[Home](#)

["Westgard Rules"](#)

[Essays](#)

[Basic QC Practices](#)

[CLIA](#)

[High Reliability](#)

["Housekeeping"](#)

[ISO](#)

[Links](#)

[Maryland General](#)

[Method Validation](#)

[Personal](#)

[QC Design](#)

[Quality Requirements and Standards](#)

[Quality of Laboratory Testing](#)

Statistics

Six Sigma

Tools

Trends

Guest Essay

Risk Management

QC Applications

QC Design

Basic QC Practices

Method Validation

Six Sigma

Sigma Metric Analysis

Quality Standards

Lessons

Basic QC Practices

Basic Planning for Quality

Basic Method Validation

Z-Stats / Basic Statistics

Quality Management

Advanced Quality Management / Six Sigma

"Westgard Rules"

Patient Safety Concepts

High Reliability

ISO

CLIA & Quality

Quality Requirements

CLIA Final Rule

[Downloads](#)

[Store](#)

[Resources](#)

[About Us](#)

[Feedback Form](#)



Westgard QC



[Home](#)

[CLIA & Quality](#)

[Quality Requirements](#)

[Desirable Biological Variation Database specifications](#)

[Quality Requirements](#)

[Desirable Biological Variation Database specifications](#)

The 2014 edition of Desirable Specifications for imprecision, inaccuracy, and total allowable error, calculated from data on within-subject and between-subject biologic variation. This database was updated and compiled by Dr. Carmen Ricos and colleagues, before the EFLM took the helm. The EFLM then instructed us to stop updating our tables. We were honored to host this database for 15 years.

Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation

This most recent and extensive listing of biologic goals has been provided by Ricos C, Alvarez V, Cava F, Garcia-Lario JV, Hernandez A, Jimenez CV, Minchinela J, Perich C, Simon M. "Current databases on biologic variation: pros, cons and progress." Scand J Clin Lab Invest 1999;59:491-500. This database was most recently updated in 2014.

Please note: Future updates of this database will be handled by the EFLM! They have revoked our permission to display these data.

Annex I, Part I: Within-subject and between-subject CV values of analytes and Desirable Analytical Quality Specifications for imprecision, bias and total error

11-Desoxycortisol through a-Fetoprotein
 Albumin through CA 549 antigen
 Calcium through Cystine
 Dehydroepiandrosterone sulfate through Hydroxylproline
 Immunoglobulin A through Lycopene
 Magnesium through Oxalate, output
 pCO₂ through Rheumatoid factor
 SCC antigen through Zinc

See The Reference List

See The References

See The original Guest Essay

Note on abbreviations:

CVI = within-subject biologic variation

CVG = between-subject biologic variation

I = desirable specification for imprecision

B = desirable specification for inaccuracy

TE = desirable specification for allowable total error

	Analyte	Number of Papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	11-Desoxycortisol	2	21.3	31.5	10.7	9.5	27.1
S-	17-Hydroxyprogesterone	2	19.6	50.4	9.8	13.5	29.7
U-	4-hydroxy-3-methoximandelate (VMA)	1	22.2	47.0	11.1	13.0	31.3
S-	5' Nucleotidase	2	23.2	19.9	11.6	7.6	26.8
U-	5'-Hydroxyindolacetate, concentration	1	20.3	33.2	10.2	9.7	26.5
S-	α1-Acid Glycoprotein	3	11.3	24.9	5.7	6.8	16.2
S-	α1-Antichymotrypsin	1	13.5	18.3	6.8	5.7	16.8

S-	α 1-Antitrypsin	3	5.9	16.3	3.0	4.3	9.2
S-	α 1-Globulins	2	11.4	22.6	5.7	6.3	15.7
U-	α 1-Microglobulin, concentration, first morning	1	33.0	58.0	16.5	16.7	43.9
P-	α 2-Antiplasmin	1	6.2	---	3.1	---	---
S-	α 2-Globulins	2	10.3	12.7	5.2	4.1	12.6
S-	α 2-Macroglobulin	4	3.4	18.7	1.7	4.75	7.56
U-	α 2-Microglobulin output, first morning	1	29.0	32.0	14.5	10.8	34.7
P-	α -aminobutyric acid	1	24.7	32.3	12.4	10.2	30.5
S-	α -Amylase	7	8.7	28.3	4.4	7.4	14.6
S-	α -Amylase (pancreatic)	2	11.7	29.9	5.9	8.0	17.7
U-	α -Amylase (pancreatic)	2	69.5	105.0	34.75	31.48	88.82
U-	α -Amylase concentration, random	1	94.0	46.0	47.0	26.2	103.7
P-	α -Carotene	1	24.0	65.0	12.0	17.3	37.1
S-	α -Carotene	1	48.0	65.0	24.0	20.2	59.8
S-	α -Fetoprotein(non hepatic carcinoma)	2	12.2	45.6	6.1	11.8	21.9
S-	α -Tocopherol	3	13.8	15.0	6.9	5.1	16.5
S-	Acid phosphatase	2	8.9	8.0	4.5	3.0	10.3
S-	Acid phosphatase tartrate-resistant (TR-ACP)	2	8.0	13.3	4.0	3.9	10.5
S-	Acid phosphatase prostatic activity (PAP)	1	33.8	---	16.9	---	---
P-	Activated partial thromboplastine time	3	2.7	8.6	1.4	2.3	4.5
P-	Adiponectin	1	18.8	51.2	9.4	13.6	29.1

S-	Adenosine deaminase (ADA)	1	11.7	25.5	5.9	7.0	16.7
P-	Alanine	1	14.7	55.8	7.4	14.4	26.6
S-	Alanine aminopeptidase	1	4.1	---	2.1	---	---
S-	Alanine aminotransferase (ALT)	9	19.40	41.6	9.7	11.48	27.48

[Back to top](#)

	Analyte	Number of papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	Albumin	24	3.2	4.75	1.6	1.43	4.07
U-	Albumin, concentration, first morning	3	36.0	55.0	18.0	16.4	46.1
U-	Albumin, output, night urine	3	29.5	58.0	14.8	16.3	40.6
S-	Albumin, glycated	1	5.2	10.3	2.6	2.9	7.2
S-	Aldosterone	2	29.4	40.1	14.7	12.4	36.7
U-	Aldosterone	1	39.4	40.1	19.7	14.05	46.56
S-	Alkaline phosphatase	22	6.45	26.1	3.23	6.72	12.04
S-	Alkaline phosphatase, bone	4	6.2	37.4	3.1	9.5	14.6
S-	Alkaline phosphatase, liver	1	10.0	27.0	5.0	7.2	15.4
S-	Alkaline phosphatase, placental	1	19.1	---	9.6	---	---
U-	Ammonia, output, 24h	1	24.7	27.3	12.4	9.2	29.6
S-	Amyloid A	1	25.0	61.0	12.5	16.5	37.1
S-	Androstendione	2	15.8	38.8	7.9	10.47	23.51
S-	Anion gap		9.5	10.1	4.8	3.5	11.3

P-	Antithrombin III	4	5.2	15.3	2.6	4.0	8.3
S-	Apolipoprotein A1	11	6.5	13.4	3.3	3.7	9.1
S-	Apolipoprotein B	9	6.9	22.8	3.5	6.0	11.6
P-	Arginine	1	19.3	34.1	9.7	9.8	25.7
S-	Arilestease activity, non inhibited	1	3.8	37.2	1.9	9.3	12.5
P-	Ascorbate (Vitamin C)	1	20.0	21.0	10.0	7.3	23.8
S-	Ascorbate (Vitamin C)	3	26.0	31.0	13.0	10.1	31.6
P-	Asparagine	1	12.3	28.0	6.2	7.6	17.8
S-	Aspartate aminotransferase (AST)	13	12.3	23.1	6.15	6.54	16.69
P-	Aspartic acid	1	31.2	55.1	15.6	15.8	41.6
S-	β -2-Microglobulin	1	5.9	15.5	3.0	4.1	9.0
P-	β -Carotene	1	18.0	48.0	9.0	12.8	27.7
S-	β -Carotene	4	36.0	39.7	18.0	13.4	43.1
S-	β -Cryptoxantin	1	36.7	---	18.4	---	---
S-	β -Globulins	2	10.1	9.1	5.1	3.4	11.7
B-	Base excess	1	76.4	43.2	38.2	21.9	85.0
B-	Basophile, count	3	28.0	54.8	14.0	15.4	38.5
S-	Bilirubin total	11	21.8	28.4	10.90	8.95	26.94
S-	Bilirubin conjugated	2	36.8	43.2	18.4	14.2	44.5
P-	C Protein	1	5.6	55.2	2.9	13.9	18.7
S-	C reactive protein	3	42.2	76.3	21.1	21.8	56.6
S-	C reactive protein - high sensitivity	1	49.70	89.23	24.85	25.53	66.54

S-	C3 Complement	2	5.2	15.6	2.6	4.1	8.4
S-	C4 Complement	2	8.9	33.4	4.5	8.6	16.0
S-	CA 125 antigen	4	24.7	54.6	12.4	15.0	35.4
S-	CA 15.3 antigen	5	6.1	62.9	3.1	15.8	20.8
S-	CA 19.9 antigen	2	16.0	130.5	7.98	32.87	46.03
S-	CA 549 antigen	1	9.1	33.4	4.6	8.7	16.2

[Back to top](#)

	Analyte	Number of papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	Calcium	24	2.1	2.5	1.05	0.82	2.55
S-	Calcium, complexed	1	5.3	4.5	2.7	1.7	6.1
U-	Calcium, concentration, 24h	4	27.5	36.6	13.8	11.4	34.1
S-	Calcium, ionized	2	1.7	1.9	0.9	0.6	2.0
S-	Calcium, protein bound	1	4.1	6.1	2.1	1.8	5.2
S-	Calcium, ultrafiltrable	1	2.2	2.7	1.1	0.9	2.7
S-	Carbohydrate deficient transferrin	1	7.1	38.7	3.6	9.8	15.7
B-	CO ₂ , total	1	4.0	4.8	2.0	1.56	4.86
S-	Carcinoembryonic antigen (CEA)	10	12.7	55.6	6.4	14.3	24.7
S-	Carnitine, Acyl-free	1	11.35	24.3	5.68	6.71	16.07
S-	Carnitine, free	1	8.05	16.65	4.03	4.62	11.26
S-	Carnitine, total	1	8.85	11.80	4.43	3.69	10.99

S-	Ceruloplasmin (ferroxidase)	2	5.8	11.1	2.9	3.1	7.9
S-	Chloride	19	1.2	1.5	0.6	0.5	1.5
S-	Cholesterol	46	5.95	15.3	2.98	4.1	9.01
S-	Cholinesterase, concentration	2	7.1	---	3.6	---	---
S-	Cholinesterase, activity	3	6.1	18.2	3.1	4.8	9.8
P-	Chromogranin A	1	12.8	26.3	6.4	7.3	17.9
P-	Citrulline	1	21.4	43.9	10.7	12.2	29.9
S-	Collagen type I C propeptide (PICP)	3	7.8	26.7	3.9	7.0	13.4
S-	Collagen type I N propeptide (PINP)	3	7.4	57.3	3.7	14.4	20.5
S-	Collagen type III N propeptide (PIIINP)	1	13.6	87.2	6.8	22.1	33.3
U-	Color, first morning	1	30.9	47.4	15.5	14.1	39.6
P-	Copper	3	8.0	19.0	4.0	5.2	11.8
S-	Copper	2	4.7	13.6	2.35	3.6	7.47
P-	Cortisol	1	21.7	46.2	10.85	12.76	30.66
S-	Cortisol	3	15.2	38.1	7.6	10.26	22.8
S-	C Peptide	3	16.6	23.2	8.3	7.1	20.8
S-	Creatine kinase (CK)	9	22.8	40.0	11.4	11.5	30.3
S-	Creatine kinase MB, %	1	6.9	48.2	3.5	10.8	16.5
S-	Creatine kinase MB, activity	4	19.7	24.3	9.9	7.8	24.1
S-	Creatine kinase MB, mass	1	18.4	61.2	9.2	14.88	30.06
S-	Creatinine	28	5.95	14.7	2.98	3.96	8.87
U-	Creatinine, concentration, 24h	8	24.0	24.5	12.0	8.6	28.4

U-	Creatinine, concentration, first morning	8	23.2	25.7	11.6	8.7	27.8
U-	Creatinine, concentration, random	8	36.3	32.4	18.2	12.2	42.1
U-	Creatinine, output, 24h	8	11.0	23.0	5.5	6.4	15.4
S-	C-Terminal telopeptide type I collagen (CTY I)	5	10.85	30.6	5.43	8.12	15.45
S-	Cyfra 21.1 Antigen	2	22.2	31.1	11.1	9.6	27.9
P-	Cystatin C	1	5.5	---	2.8	---	---
S-	Cystatin C	4	5.0	13.0	2.5	3.48	7.61
P-	Cysteine	1	5.9	12.3	3.0	3.4	8.3
P-	Cystine	1	38.3	48.5	19.2	15.4	47.0

[Back to top](#)

	Analyte	Number of Papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	D-Dimer (MoM)	1	23.3	26.5	11.65	8.82	28.04
S-	Dehydroepiandrosterone sulfate (DHEAS)	3	6.35	30.70	3.188	7.84	13.08
U-	Deoxypyridinoline/creatinine, 24h	2	16.0	30.7	8.0	8.7	21.9
U-	Deoxypyridinoline/creatinine, first morning	2	13.8	34.6	6.9	9.3	20.7
U-	Deoxypyridinoline/minute, first morning	2	15.4	30.3	7.7	8.5	21.2
P-	Dipeptidyl-peptidase IV (ACE)	2	8.2	14.5	4.1	4.2	10.9
S-	Dipeptidyl-peptidase IV (ACE)	1	12.5	27.7	6.3	7.6	17.9

P-	Elastase	1	12.4	15.1	6.2	4.88	15.11
B-	Eosinophils, count	3	21.0	76.4	10.5	19.8	37.1
(B)Plat-	Epinephrine	1	25.3	---	12.7	---	---
P-	Epinephrine	1	48.3	---	24.2	---	---
B-	Erythrocytes, count	7	3.2	6.3	1.6	1.7	4.4
B-	Erythrocyte distribution wide		3.5	5.7	1.8	1.7	4.6
U-	Estradiol	1	30.4	---	15.2	---	---
S-	Estradiol	5	22.5	24.4	11.25	8.3	26.86
S-	Estradiol, free	1	22.8	---	11.40	---	---
U-	Estradiol, free	1	38.6	---	19.3	---	---
P-	Factor V coagulation	1	3.6	---	1.8	---	---
P-	Factor VII coagulation	2	6.8	19.4	3.4	5.1	10.7
P-	Factor VIII coagulation	2	4.8	19.1	2.4	4.9	8.9
P-	Factor X coagulation	1	5.9	---	3.0	---	---
S-	Ferritin	6	14.2	15.0	7.1	5.2	16.9
P-	Fibrinogen	5	10.7	15.8	5.4	4.8	13.6
(B)Erthry-	Folate	1	12.0	66.0	6.0	16.8	26.7
S-	Folate	1	24.0	73.0	12.0	19.2	39.0
S-	Follicle stimulating hormone (FSH)	5	11.0	47.2	5.5	12.12	21.19
S-	Fructosamine	3	3.4	5.9	1.7	1.7	4.5
S-	Galactosyl hydroxylysine	1	11.8	25.8	5.9	7.1	16.8
P-	γ-Fibrinogen	1	14.1	27.25	7.05	7.67	19.3

S-	γ-Globulins	2	14.6	12.3	7.3	4.8	16.8
S-	γ-glutamyltransferase (GGT)	10	13.4	42.15	6.7	11.06	22.11
S-	Globulins, total	1	5.5	12.9	2.8	3.5	8.0
P-	Glucose	1	4.5	5.8	2.3	1.8	5.5
S-	Glucose	15	5.6	7.5	2.8	2.34	6.96
(B)Erythr-	Glucose-6-phosphate-1-dehydrogenase (G6PDH)	1	32.8	31.8	16.4	11.4	38.5
B - spot	Glucose-6-phosphate-1-dehydrogenase (G6PDH)	1	7.3	10.3	3.7	3.2	9.2
P-	Glutamic acid	1	46.4	79.9	23.2	23.1	61.4
P-	Glutamine	1	12.1	22.0	6.1	6.3	16.3
S-	Glutathion peroxidase	1	7.2	21.7	3.6	5.7	11.7
P-	Glycine	1	11.8	40.3	5.9	10.5	20.2
P-	Haptoglobin	1	20.0	27.9	10.0	8.6	25.1
S-	Haptoglobin	3	20.4	36.4	10.2	10.4	27.3
S-	HDL cholesterol	25	7.3	21.2	3.65	5.61	11.63
S-	HDL 1 cholesterol	1	5.5	27.2	2.8	6.9	11.5
S-	HDL 2 cholesterol	6	15.7	40.7	7.9	10.9	23.9
S-	HDL 3 cholesterol	6	7.0	14.3	3.5	4.0	9.8
B-	Hematocrit	11	2.7	6.41	1.35	1.74	3.97
B-	Hemoglobin	13	2.85	6.8	1.43	1.84	4.19
B-	Hemoglobin A1 C	8	1.9	5.7	0.9	1.5	3.0
B-	Hemoglobin A2	1	0.7	7.7	0.35	1.93	2.51

P-	Histidine	1	9.7	27.2	4.9	7.2	15.2
P-	Homocysteine	3	8.3	33.5	4.15	8.63	15.48
S-	Hyaluronic acid	1	62.00	---	31.00	---	---
S-	Hydroxybutyrate dehydrogenase	1	6.6	---	3.3	---	---
P-	Hydroxyproline	1	34.5	56.7	17.3	16.6	45.1
U-	Hydroxyproline/minute, first morning	1	36.1	38.8	18.1	13.2	43.0
U-	Hydroxyproline/minute, second void	1	40.5	32.9	20.3	13.0	46.5

[Back to top](#)

	Analyte	Number of papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	Immunoglobulin A	7	5.4	35.9	2.7	9.1	13.5
S-	Immunoglobulin G	6	4.5	16.5	2.3	4.3	8.0
S-	Immunoglobulin M	5	5.9	47.3	3.0	11.9	16.8
S-	Immunoglobulins κ chains	3	4.8	15.3	2.4	4.0	8.0
S-	Immunoglobulins λ chains	3	4.8	18.0	2.4	4.7	8.6
S-	Inhibin B	1	10.0	25.0	5	6.73	14.98
S-	Insulin	4	21.1	58.3	10.6	15.5	32.9
S-	Insulin-like growth factor (IGF-1)	2	14.6	45.4	7.3	11.9	24.0
S-	Insulin-like growth factor binding protein 3 (IGFBP-3)	1	10.1	63.9	5.1	16.2	24.5
S-	Intercellular adhesion molecule-1 (ICAM-1)	1	1.9	21.0	1.0	5.3	6.8

(B)Leuc-	Interferon receptor	1	14.0	20.0	7.0	6.1	17.7
S-	Interleukin 1- β	1	30.0	36.0	15.0	11.7	36.5
S-	Interleukin-8	1	24.0	31.0	12.0	9.8	29.6
S-	Iron	11	26.5	23.2	13.3	8.8	30.7
P-	Isoleucine	1	15.5	45.5	7.8	12.0	24.8
S-	Kallicrein 6	1	11.80	27.6	5.9	7.5	17.24
B-	Lactate	1	27.2	16.7	13.6	8.0	30.4
S-	Lactate dehydrogenase (LDH)	11	8.6	14.7	4.3	4.3	11.4
S-	Lactate dehydrogenase 1 isoform (LDH1)	2	2.3	8.3	1.2	2.2	4.1
S-	Lactate dehydrogenase 2 isoform (LDH2)	1	3.3	2.4	1.7	1.0	3.7
S-	Lactate dehydrogenase 3 isoform (LDH3)	1	2.8	3.8	1.4	1.2	3.5
S-	Lactate dehydrogenase 4 isoform (LDH4)	1	5.9	5.3	3.0	2.0	6.9
S-	Lactate dehydrogenase 5 isoform (LDH5)	1	8.0	9.6	4.0	3.1	9.7
P-	Lactoferrin	1	11.8	23.7	5.9	6.6	16.4
S-	LDL Cholesterol	6	7.8	20.4	3.9	5.46	11.9
P-	LDL Cholesterol (oxidized)	1	21.0	50.0	10.5	13.6	30.9
S-	LDL Cholesterol, small dense	1	9.1	20.0	4.55	5.49	13.0
S-	LDL receptor mRNA	1	21.5	13.6	10.8	6.4	24.1
P-	Leucine	1	14.8	44.0	7.4	11.6	23.8

B-	Leukocytes count	8	11.4	21.3	5.73	6.05	15.49
S-	Lipase	3	32.2	31.8	16.1	11.31	37.88
S-	Lipoprotein (a)	3	20.8	18.1	10.4	6.9	24.1
P-	Lutein	1	13.0	21.0	6.5	6.2	16.9
S-	Lutein	1	23.7	---	11.9	---	---
S-	Luteinizing hormone (LH)	5	23.0	27.4	11.5	8.94	27.92
P-	Lycopene	1	22.0	33.0	11.0	9.9	28.1
S-	Lycopene	2	40.1	33.0	20.1	13.0	---
B-	Lymphocytes, count	1	10.2	35.3	5.1	9.19	17.6
B-	Lymphocytes CD4	5	25.0	---	12.5	---	---
P-	Lysine	1	11.5	38.2	5.8	10.0	19.5

[Back to top](#)

	Analyte	Number of Papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
(B)Erythr-	Magnesium	2	5.6	11.3	2.8	3.2	7.8
(B)Leuc-	Magnesium	1	18.5	12.4	9.25	5.57	20.83
(B)Mon -	Magnesium	1	18.1	20.3	9.1	6.8	21.7
S-	Magnesium	9	3.6	6.4	1.8	1.8	4.8
U-	Magnesium, ionized	1	1.9	5.1	1.0	1.4	2.9
U-	Magnesium, output, 24h	2	38.3	37.6	19.2	13.4	45.0
(B)Erythr-	Mean corpuscular hemoglobin (MCH)	4	1.4	5.2	0.7	1.35	2.5

(B)Erythr-	Mean corpuscular hemoglobin concentration (MCHC)	5	1.06	1.2	0.53	0.4	1.27
(B)Erythr-	Mean corpuscular volume (MCV)	7	1.4	4.85	0.7	1.26	2.42
(B)Plat-	Mean platelet volume (MPV)	3	4.3	8.1	2.15	2.29	5.84
P-	Metionine	1	14.7	43.4	7.4	11.5	23.6
B-	Monocytes, count	3	17.8	49.8	8.9	13.2	27.9
S-	Myeloperoxidase	1	36.0	30.0	18.0	11.7	41.4
S-	Myoglobin	2	13.9	29.6	7.0	8.2	19.6
U-	N-Acetyl Glucosaminidase, concentration, first morning	2	52.9	22.0	26.5	14.3	58.0
U-	N-Acetyl Glucosaminidase/Creatinine	2	51.1	21.8	25.6	13.9	56.0
B-	Neutrophyles, count	5	17.1	32.8	8.55	9.25	23.35
U-	Nitrogen, output	1	13.9	24.2	7.0	7.0	18.4
B(Plat)-	Norepinephrine	1	9.5	---	4.8	---	---
P-	Norepinephrine	1	19.5	---	9.8	---	---
U-	N-Telopeptide type I collagen	3	15.5	37.6	7.75	10.17	22.95
S-	N-terminal (NT)-proBNP	2	10.0	16.0	5.0	4.7	13.0
P-	Ornithine	1	18.4	54.9	9.2	14.5	29.7
P-	Osmolality	1	1.3	1.5	0.7	0.5	1.6
Saliva-	Osmolality	1	9.5	35.8	4.8	9.3	17.1
S-	Osmolality	1	1.3	1.2	0.7	0.4	1.5
U-	Osmolality, first morning	1	28.3	57.9	14.2	16.1	39.5
S-	Osteocalcin	5	6.35	30.9	3.18	7.89	13.13

U-	Oxalate, concentration, 24h	1	44.0	18.0	22.0	11.9	48.2
U-	Oxalate, output, 24h	1	42.5	19.9	21.3	11.7	46.8

Back to top

	Analyte	Number of papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
B-	pCO ₂	1	4.8	5.3	2.4	1.8	5.7
B-	pH [H ⁺]	1	3.5	2.0	1.8	1.0	3.9
B-	pH (pH units)	1	0.2	---	0.1	---	---
S-	Paraoxonase 1	1	13.4	84.0	6.7	21.3	32.3
S-	Paraoxonase 1 substrate inhibition (PON 4SI)	1	3.9	80.1	1.9	20.0	23.2
S-	Paraoxonase, activity (salt stimulated)	1	8.0	86.4	4.0	21.7	28.3
P-	Parathyroid hormone (PTH)	1	25.3	43.4	12.65	12.56	33.43
S-	Parathyroid hormone (PTH)	1	25.9	23.8	13.0	8.8	30.2
S-	Phenylacetate	1	6.6	25.2	3.3	6.5	12.0
P-	Phenylalanine	1	9.5	40.6	4.8	10.4	18.3
S-	Phosphate	17	8.15	10.8	4.08	3.38	10.11
U-	Phosphate, output, 24h	5	18.0	22.6	9.0	7.2	22.1
Patient-	Phosphate tubular reabsorption	1	2.7	3.3	1.4	1.1	3.3
S-	Phospholipids	1	6.5	11.1	3.3	3.2	8.6
P-	Plasminogen	1	7.7	---	3.9	---	---
B-	Platelets, count	7	9.1	21.9	4.6	5.9	13.4

B-	Platelet distribution wide	2	2.8	---	1.4	---	---
B-	Plateletcrit	2	11.9	---	6.0	---	---
U-	Porphobilinogen	1	17.0	31.0	8.5	8.8	22.9
U-	Porphyrins (total)	1	40.0	---	20.0	---	---
(B)Leuc-	Potassium	1	13.6	13.4	6.8	4.8	16.0
S-	Potassium	20	4.6	5.6	2.3	1.81	5.61
U-	Potassium, output	4	24.4	22.2	12.2	8.2	28.4
S-	Prealbumin	1	10.9	19.1	5.5	5.5	14.5
S-	Pregnancy-associated plasma protein A (PAPP-A)	1	12.6	14.0	6.3	4.71	15.1
P-	Prolactin	1	39.2	65.1	19.6	19.0	51.34
S-	Prolactin	4	23.0	35.0	11.5	10.5	29.4
P-	Proline	1	17.0	104.4	8.5	26.4	40.5
P-	Prolyl endopeptidase	2	16.8	13.9	8.4	5.5	19.3
S-	Properdin factor B	1	9.5	11.2	4.7	3.7	11.5
S-	Prostatic specific antigen (PSA)	3	18.1	72.4	9.1	18.7	33.6
S-	Protein	18	2.75	4.7	1.38	1.36	3.63
S-	Protein, glycated	1	0.9	11.6	0.5	2.9	3.7
U-	Protein, output, 24h	2	35.5	23.7	17.8	10.7	40.0
P-	Prothrombin time	2	4.0	6.8	2.0	2.0	5.3
U-	Pyridinoline	1	19.4	23.6	9.7	7.6	23.6
B-	Pyruvate	1	15.2	13.0	7.6	5.0	17.5

S-	Receptor for advanced glycation end-products (RAGE)	1	14.6	56.5	7.3	14.59	26.63
B-	Red cell distribution wide (RDW)	4	3.5	5.7	1.8	1.7	4.6
B-	Reticulocyte highly fluorescent, count	1	10.0	62.0	5.0	15.7	24.0
B-	Reticulocyte low fluorescent, count	1	1.6	4.9	0.8	1.3	2.6
B-	Reticulocyte medium fluorescent, count	1	13.0	33.0	6.5	8.9	19.6
B-	Reticulocyte, count	1	11.0	29.0	5.5	7.8	16.8
P-	Retinol	1	6.2	21.0	3.1	5.5	10.6
S-	Retinol	2	13.6	19.0	6.8	5.8	17.1
S-	Rheumatoid factor	1	8.5	24.5	4.3	6.5	13.5

[Back to top](#)

<

	Analyte	Number of papers	Biological Variation		Desirable specification		
			CVI	CVg	I(%)	B(%)	TE(%)
S-	SCC antigen	1	39.4	35.7	19.7	13.3	45.8
P-	S Protein	1	5.8	63.4	2.9	15.9	20.7
P-	Selenium	1	12.0	14.0	6.0	4.6	14.5
B-	Selenium	1	12.0	12.0	6.0	4.2	14.1
P-	Serine	1	12.8	42.8	6.4	11.2	21.7
S-	Sex hormone binding globulin (SHBG)	2	13.05	36.35	6.53	9.66	20.42
(B)Erythr-	Sodium	1	1.8	12.4	0.9	3.1	4.6
(B)Leuc-	Sodium	1	51.0	36.4	25.5	15.7	57.7

S-	Sodium	21	0.6	0.7	0.3	0.23	0.73
B-	Sodium Bicarbonate	1	4.0	4.8	2.0	1.6	4.9
S-	Sodium Bicarbonate	7	4.8	4.7	2.4	1.7	5.6
Sweat-	Sodium Chloride	1	15.0	25.0	7.5	7.3	19.7
U-	Sodium, output, 24 h.	4	28.7	16.7	14.4	8.3	32.0
P-	Soluble CD163	1	9.0	35.9	4.5	9.3	16.7
U-	Specific gravity	1	0.4	1.0	0.2	0.27	0.60
Semen-	Spermatozoa, concentration	1	26.8	56.4	13.4	15.6	37.7
Semen-	Spermatozoa, morphology	1	19.6	44.0	9.8	12.0	28.2
Semen-	Spermatozoa, progressive motility	1	15.2	32.8	7.6	9.0	21.6
Semen-	Spermatozoa, fast progressive motility	1	18.8	51.8	9.4	13.8	29.3
Semen-	Spermatozoa, total motility	1	18.4	29.8	9.2	8.8	23.9
Semen-	Spermatozoa, vitality	1	10.3	25.8	5.2	6.9	15.4
S-	Superoxide dismutase	1	17.1	10.5	8.6	5.0	19.1
(B)Erythr-	Superoxide dismutase	1	12.3	4.9	6.2	3.3	13.5
P-	Taurine	1	30.6	44.0	15.3	13.4	38.6
P-	Testosterone	1	12.6	40.80	6.3	10.68	21.07
S-	Testosterone	7	9.25	22.05	4.63	5.98	13.61
Saliva-	Testosterone	1	17.3	28.8	8.7	8.4	22.7
U-	Testosterone	1	25.0	---	12.5	---	---
S-	Testosterone, free	3	9.3	---	4.7	---	---
U-	Testosterone, free	1	51.7	---	25.9	---	---

S-	Thyroglobulin	2	14.0	39.0	7.0	10.4	21.9
S-	Thyroglobulin antibody	1	8.5	82.0	4.3	20.6	27.6
S-	Thyroid peroxidase antibody	1	11.3	147.0	5.7	36.9	46.2
P-	Thyroid stimulating hormone (TSH)	1	29.30	48.4	14.65	14.14	38.2
S-	Thyroid stimulating hormone (TSH)	9	19.3	24.6	9.7	7.8	23.7
S-	Thyrotropin receptor antibody	1	4.8	---	2.4	---	---
S-	Thyroxine (T4)	11	4.9	10.9	2.5	3.0	7.0
P-	Thyroxine, free (FT4)	1	7.1	9.1	3.55	2.89	8.74
S-	Thyroxine, free (FT4)	5	5.7	12.1	2.9	3.3	8.0
S-	Thyroxine binding globulin (TBG)	2	0.09	0.06	0.0	0.0	0.1
P-	Tirosine	1	10.5	61.0	5.3	15.5	24.1
S-	Tissue polypeptide antigen (TPA)	1	31.1	63.7	15.6	17.7	43.4
S-	Tissue polypeptide specific antigen (TPS)	1	36.1	108.0	18.1	28.5	58.3
U-	Total catecholamines, concentration, 24h	1	24.0	32.0	12.0	10.0	29.8
S-	Transferrin	5	3.0	4.3	1.5	1.3	3.8
P-	Treonine	1	17.9	33.1	9.0	9.4	24.2
S-	Triglyceride	31	19.9	32.7	9.95	9.57	25.99
P-	Triiodothyronine (T3)	1	9.4	18.5	4.7	5.19	12.94
S-	Triiodothyronine (T3)	10	6.9	12.3	3.45	3.53	9.22
S-	Triiodothyronine, free (FT3)	4	7.9	17.6	4.0	4.8	11.3
P-	Troponin I	1	37.1	179.2	18.55	45.75	76.36

S-	Troponin I	5	14.05	63.75	7.03	16.32	27.91
S-	Troponin T	1	30.5	90.0	15.3	23.7	48.9
P-	Tryptophan	1	22.7	152.6	11.4	38.6	57.3
S-	Tumor Necrosis Factor- α (TNF- α)	1	43.0	29.0	21.5	13.0	48.4
S-	Urate	16	8.6	17.5	4.3	4.87	11.97
U-	Urate, output, 24h	2	16.8	14.4	8.4	5.53	19.39
S-	Urea	20	12.1	18.7	6.05	5.57	15.55
U-	Urea, output, 24h	4	17.4	25.4	8.7	7.7	22.1
P-	Valine	1	10.6	40.1	5.3	10.4	19.1
U-	Vanilmandelic Acid concentration, 24h	1	22.2	47.0	11.1	13.0	31.3
S-	Vascular cell adhesion molecule-1 (VCAM-1)	1	5.2	16.0	2.6	4.2	8.5
P-	Vascular endothelial growth factor	1	14.1	18.1	7.1	5.7	17.4
B-	Vascular endothelial growth factor	1	14.3	28.8	7.2	8.0	19.8
S-	Vascular endothelial growth factor	1	10.7	47.6	5.4	12.2	21.0
P-	Vitamin B1	1	4.8	12.0	2.4	3.2	7.2
B-	Vitamin B2 (Riboflavin)	1	5.8	10.0	2.9	2.9	7.7
(B)Eryth-	Vitamin B2 (Riboflavin)	1	6.4	11.0	3.2	3.2	8.5
(B)Eryth-	Vitamin B2 status (gluthation reductase activation)	1	5.2	40.0	2.6	10.1	14.4
(B)Eryth-	Vitamin B12	1	15.0	69.0	7.5	17.7	30.0
(B)Eryth-	Vitamin B6	1	14.0	24.0	7.0	6.9	18.5
B-	Vitamin B6	1	20.0	34.0	10.0	9.9	26.4

(B)Eryth-	Vitamin E (Tocopherol)	1	7.6	21.0	3.8	5.6	11.9
(B)Eryth-	Vitamin K (Phylloquinone)	1	38.0	44.0	19.0	14.5	45.9
S-	VLDL Cholesterol	2	27.6	---	13.8	---	---
P-	Von Willebrand factor	3	2.5	27.3	1.3	6.9	8.9
S-	Water	1	3.1	0.1	1.6	0.8	3.3
S-	Zeaxanthine	1	34.7	---	17.4	---	---
S-	Zinc	1	9.3	9.4	4.7	3.3	11.0
P-	Zinc	3	11.0	14.0	5.5	4.5	13.5

[Back to top](#)

[Joomla SEF URLs by Artio](#)



JAMES WESTGARD
FOUNDER

[Blog](#)

[About Us](#)

[Reference Materials
& Resources](#)

[Calculators](#)

[QC Tools](#)

[QC Calculators](#)

[Method Validation Tools](#)

[Six Sigma Calculators](#)

[Normalized OPSpecs Calculator](#)

[Quality Control Grid Calculator](#)

[Control Limit Calculator](#)

[Reportable Range Calculator: Quantifying Errors](#)

[Reportable Range Calculator: Recording Results](#)

[Dispersion Calculator and Critical Number of Test Samples](#)

[Online Store](#)



[Nothing but the Truth Manual](#)

\$225.00



[Six Sigma Risk Analysis](#)

\$90.00

[Photo Gallery](#)



WHAT'S POPULAR

WHAT'S NEW

Member Login

Øverst i skjemaet

To access the private area of this site, please log in.

Remember Me ☐

[Log in](#)

[Forgot Login?](#)

[Sign up](#)

[Nederst i skjemaet](#)

[Other Articles](#)

[1999 Stockholm Consensus Statement](#)

[2016 State of the Art Hematology Performance Specifications](#)

[A selection of SEKK-DMax specifications](#)

[Belgium EQA performance specifications](#)

[Biological Variation Database reference list](#)

[Biological Variation Database references](#)

[Biological Variation in Patients with Disease](#)

[CLIA Requirements for Analytical Quality](#)

[Clinical Quality Requirements](#)

[Consolidated Comparison of Hematology Performance Specifications](#)

[Desirable Biological Variation Database specifications](#)

[European Biologic Goals](#)

[Final Milan Hierarchy](#)

[French EQA performance specifications](#)

[Milan 2014 Consensus Draft on Quality Specifications](#)

[Minimum analytical quality requirements](#)

[Minimum Specifications from Biological Variation database](#)

[Optimal Biological Variation database specifications](#)

[Quality Requirements for Dogs, Cats, and Horses?](#)

[RCPA Allowable Limits of Performance for Biochemistry](#)

[Rilibak - German Guidelines for Quality](#)

[Spanish Minimum Consensus Performance Specifications](#)

[The New Era of Biological Variation: EUBIVAS](#)

[Øverst i skjemaet](#)

[Nederst i skjemaet](#)

[What's New](#)

[West China Hospital Sigma Verification](#)

[The 17 Signs of a Struggling \(Losing\) Laboratory](#)

[Affiliated Hospital of Qingdao University - Sigma Verification of Performance](#)

[Cookie Policy for Westgard QC](#)

[Privacy Policy](#)



Copyright © 2009. All rights reserved.

Westgard QC • 7614 Gray Fox Trail • Madison, Wisconsin 53717

Call 608-833-4718 or E-mail westgard@westgard.com

Joomla Web Design by IDL Web Inc Joomla Hosting by IDL Web Host

"Westgard Rules"

[Questions](#)

[Interviews](#)

[Lessons](#)

[CLIA & Quality](#)

[Essays](#)

[Tools](#)

[QC Applications](#)

[Photos](#)

[Contact WQC](#)

[Site Map](#)

[Home](#)

["Westgard Rules"](#)

[Essays](#)

[Basic QC Practices](#)

[CLIA](#)

[High Reliability](#)

["Housekeeping"](#)

[ISO](#)

[Links](#)

[Maryland General](#)

[Method Validation](#)

[Personal](#)

[QC Design](#)

[Quality Requirements and Standards](#)

[Quality of Laboratory Testing](#)

[Statistics](#)

[Six Sigma](#)

Tools

Trends

Guest Essay

Risk Management

QC Applications

QC Design

Basic QC Practices

Method Validation

Six Sigma

Sigma Metric Analysis

Quality Standards

Lessons

Basic QC Practices

Basic Planning for Quality

Basic Method Validation

Z-Stats / Basic Statistics

Quality Management

Advanced Quality Management / Six Sigma

"Westgard Rules"

Patient Safety Concepts

High Reliability

ISO

CLIA & Quality

Quality Requirements

CLIA Final Rule

Downloads

Store

[Resources](#)

[About Us](#)

[Feedback Form](#)

[loading](#)